



ZENITH HUMAN MILK BANK EQUIPMENTS

**An ISO 9001:2015, ISO 13485:2016,
CE & GMP Certified**

Website: www.sewardsystems.com

Milk Containers

BPA Free Polypropylene Bottle



Available single packed or bulk with cap, option to include test.

- Sizes: 50ml, 130ml, 240ml
- Sterile or gamma disinfection options

Features & Benefits

- Suitable for the collection, storage and pasteurisation of breast milk or formula
- Packaged sterile and ready to use or gamma disinfected
- Freezable down to -20°C
- EU Standard bottle thread
- Calibrated graduation (meets the requirements of the standard EN 14350-1)
- Suitable for standard bottle warmers and transport baskets
- Leak proof due to seamless manufacturing

ZENITH BPA Free Milk Containers



Specifications :

- Material Used SS 304
- BPA Free Milk Containers
- Milk Storage Capacity : 200ml/250ml.

Bottle Sealer with Foil for sealing Bottles



Specifications:

1. Specifically designed to avoid leakage of tank water into the breast milk bottle during the heating and cooling cycles of the breast milk
2. This equipment is intended to seal metal foiled wafer inside the capped containers. Heating takes place in the metal foil and conducts heat to its plastic coating and subsequently causes the.
3. The package should consist of the induction heating generator
4. Aluminium Foil Sealer
5. Container Sealing Machines
6. Power : 200 W
7. Foil Sealing : 2000nos (free)
8. Sealing width : 175mm
9. With digital display.



Save Milk Save Life...!

Marketed By: Seward Systems Inc, Maharashtra

India - 400 066,

Service support Email ID: info@sewardsystems.com

Mobile. No +91- 9820426021, 9820426011