



ZENITH HUMAN MILK BANK EQUIPMENTS

An ISO 9001:2015, ISO 13485:2016,
CE & GMP Certified

Website: www.sewardsystems.com

ZENITH Deep Freezer & Refrigerator

DF 325R & RF 325

This Special Deep Freezer & Refrigerators are only made For Human Milk Banking. It's made from fully food grade Stainless steel materials. This special Deep freezer is made with Rapid cooling technology



Optional features:

- Food Grade Fully SS Construction.
- Vertical deep freezer with double door glass having temperature -18°C to -24°C .
- Vertical Refrigerator with double door glass having temperature 2°C to 8°C .
- Double door with Vacuum glass.
- Total 4 chambers.
- 200 bottles storage capacity.
- Fast cooling Technology.
- 300 litres capacity.
- Data Logger system (for Deep Freezer only) with pen drive for data storing & tracking.
- It has microprocessor control for operation with digital monitoring display.
- Audio & Visual high and Low temperature alarms
- Thermostat with digital Display.
- It is lockable and the door walls has insulated with sufficient thickness of high quality PUF insulation.
- CFC free refrigerants.
- 230v, 50Hz single phase AC supply.
- 75 mm puff insulation CFC frees R141B as blowing agent, density 30/35 Kg/Cu. m.
- Two models available Glass doors and without Glass full SS doors.
- ISO 9001:2015, ISO 13485:2016, CE, GMP RoHS Certification.



Save Milk Save Life...!

Marketed By: Seward Systems Inc, Maharashtra

India - 400 066,

Service support Email ID: info@sewardsystems.com

Mobile. No +91- 9820426021, 9820426011