



ZENITH HUMAN MILK BANK EQUIPMENTS

**An ISO 9001:2015, ISO 13485:2016,
CE & GMP Certified**

Website: www.sewardsystems.com

ZENITH Milk Warmer

MW 37



Specifications:

- Constructed in fully SS.
- Digital PID Controller
- Fully automatic with Timer
- 4 Bottles tray (150ml/250bottle)
- Product sensor for temperature
- Power supply 230v/50Hz

Product dimensions:

- Outer body: l-600mm x d -550mm x h – 650mm
- INNER TANK: L- 400MM X D- 350MM X H – 300MM

ZENITH Autoclave

Autoclave ACL 25



Specifications:

- These double walled units have inner chamber (Boiler) made of SS.
- Outer wall is covered with stainless steel sheet.
- Thick stainless-steel lid is tightened by radial locking system and can be lifted through a pedal lifting device for user comfort
- Fitted with joint less silicon rubber gasket.
- Heating is done through specially designed water heaters of suitable wattage to ensure optimum sterilization time in every sterilization cycle.
- Low water device cuts off electric supply to heater in case of low water level detection.
- Dimensions: 12" Dia X 22" Height.
- Power Source: 220/240 volts A/C, 50 Hz Electrical Supply, Working Pressure: 15 to 20 psi,
- Working Temperature: 121° C to 125° C, Controller: Auto set 15 psi to 121° C.
- Ideal Use: provides a physical method of sterilization by killing bacteria, viruses, and even spores present in the material put inside of the vessel using steam under pressure. used in hospitals,
- Bottle Container – 1nos. (with each autoclave)



Save Milk Save Life...!

Marketed By: Seward Systems Inc, Maharashtra

India - 400 066,

Service support Email ID: info@sewardsystems.com

Mobile. No +91- 9820426021, 9820426011