

Catheter Mounts


Catheter Mounts

Altech range of catheter mounts includes a wide choice of smoothbore PVC and extendible configurations

- Smoothbore products provide low resistance to flow while external coil reduces kinks and blockages. Besides they prevent fluid build up reducing risks of bacterial growth thanks to smooth inner surface.
- Extendible PP product provides a lightweight and flexible solution and allows clinician to adjust the position of the catheter mount
- 360 degree double swivel connector increases ease of use
- Double cap helps to maintain PEEP during suctioning or bronchoscopy.
- Luer lock port for gas sampling
- 15M or 22F adapters available




CATHETER MOUNT EXTENDIBLE



22 mmF

Dead space (ml): Close: 25, Open: 40



15 mmM


22 mmM/15mmF

Dead space (ml): Close: 25, Open: 40

Extendible


Length	Code	Qty/Box
15 cm/6"	AL-17300	50
15 cm/6"	AL-13300	50

CATHETER MOUNT SMOOTH BORE



15 mmF

Dead space (ml): 30



22 mmF

22 mmM/15mmF

Dead space (ml): 30


Smoothbore

Length	Code	Qty/Box
15 cm/6"	AL-12300	50
15 cm/6" (PVC-free)	AL-P12300	50

Smoothbore

Length	Code	Qty/Box
15 cm/6"	AL-16300	50
15 cm/6" (PVC-free)	AL-P16300	50

CATHETER MOUNT & FILTER



22 mmF

22 mmM/15mmF

Dead space (ml): Close: 103, Open: 117


Extendible

Length	Code	Qty/Box
15 cm/6"	AL-17380	20


Catheter Mounts




CATHETER MOUNTS




22 mmF




15 mmM/22 mmF



22 mmF



15 mmM/22 mmF



EXTENDIBLE TUBE

CODE: 67300
QTY/Box: 50
Dead space (ml):
Close: 22, Open: 37

SMOOTHBORE TUBE

CODE: 66300
PVC-FREE CODE: P66300
QTY/Box: 50
Dead space (ml): 27

SWIVEL CATHETER MOUNT

CODE: AL-10000
QTY/Box: 100
Dead space (ml): 3


CATHETER MOUNT SMOOTHBORE

Smoothbore


Length	Code	Qty/Box
15 cm/6"	AL-22300	50
15 cm/6" (PVC-free)	AL-P22300	50

Smoothbore

Length	Code	Qty/Box
15 cm/6"	AL-26300	50
15 cm/6" (PVC-free)	AL-P26300	50



22 mmF



22 mmM/15 mmF

22 mmF


22 mmM/15 mmF

Dead space (ml): 30

CATHETER MOUNT EXTENDIBLE

Extendible

Length	Code	Qty/Box
15 cm/6"	AL-27300	50



22 mmF

22 mmM/15 mmF

Dead space (ml): Close: 25, Open: 40