



SCEPTRE

your kind of care



Sceptre Medical Devices Pvt Ltd

CORP OFFICE : 06-Central Market, Ashok Vihar, PH-I, Delhi-110 052 India TEL : +91 11 27303001
WORKS : 187, Sector-38, PH-I, HSIIDC, Rai, Sonapat-131 029, HR, India TEL : +91 130 2368187
info@sceptremedica.com | www.sceptremedica.com

SCEPTRE

REIN
CATH


Long Term Hemodialysis Catheters

"Catheters and key components are developed based on research & are sourced from TAIWAN & PRC"

Clinically Tested

ISO 9001 : 2015 Certified

FDA Registered



Soft Tip

Reduces trauma to vessel wall upon introduction and in-situ.

Antimicrobial Heparin Cap

Specially-made adaptor helps to avoid infection.

Rotating Suture Wing

Rotating Suture Wing provides secure external anchoring.

Radiopaque Material

Sufficiently Radiopaque materials ensures the correct placement of catheter when visualized under X-Ray.

Biocompatible Polyurethane

Material tested for toxicity and hemocompatibility which becomes softer at body temperature.

Long-Term Double Lumen Catheter

Product Code	Size	Standard Packing Box/Carton	Dimensions Box/Carton
D301.0317S	12.5 Fr x 28.0 cm	05/30	318x235x128mm/490x330x400mm
D301.0327S	12.5 Fr x 32.0 cm	05/30	318x235x128mm/490x330x400mm
D301.0337S	12.5 Fr x 36.0 cm	05/30	318x235x128mm/490x330x400mm
D301.0347S	12.5 Fr x 40.0 cm	05/30	318x235x128mm/490x330x400mm
D301.0329S	14.5 Fr x 32.0 cm	05/30	318x235x128mm/490x330x400mm
D301.0339S	14.5 Fr x 36.0 cm	05/30	318x235x128mm/490x330x400mm
D301.0349S	14.5 Fr x 40.0 cm	05/30	318x235x128mm/490x330x400mm

	Tip to Hub	Tip to Cuff	End to End
25cm (Tip to Cuff)	24	19	36
28cm (Tip to Hub)	28	23	40
	32	27	44

CONTENTS

One Polyurethane Radiopaque Cuffed Catheter with Injection Site Caps
One Nitinol Core Guide Wire With Advancer

Two Dilator
One No. 11 Scalpel Blade
One Introducer Needle
One 10 ml Syringe
Two Heparin Locks

One IV Dressing
One Surgical Drape
One Mersilk Suture
One Sheath/Sheath Valved
One Tunneler (Trocar)
Four Gauze

SCEPTRE

REIN
CATH

Double & Triple Lumen Hemodialysis Catheters

"Catheters and key components are developed based on research & are sourced from TAIWAN & PRC"

Clinically Tested

ISO 9001 : 2015 Certified

FDA Registered

Soft Tip

Reduces trauma to vessel wall upon introduction and in-situ.

Antimicrobial Heparin Cap

Specially-made adaptor helps to avoid infection.

Rotating Suture Wing

Rotating Suture Wing provides secure external anchoring.

Radiopaque Material

Sufficiently Radiopaque materials ensures the correct placement of catheter when visualized under X-Ray.

Biocompatible Polyurethane

Material tested for toxicity and hemocompatibility which becomes softer at body temperature.

Double Lumen Catheter Subclavin & Curved Extensions

PRODUCT CODE					
Subclavin Extensions	Curved Extensions	Size	Standard Packing Box/Carton	Dimensions Box/Carton	Gross Wt.
D301.0351S	D301.0351C	6.5 Fr x 8-10 cm	05/30	318x235x128mm/490x330x400mm	5.0 Kg
D301.0312S	D301.0312C	8.5 Fr x 10-11 cm	05/30	318x235x128mm/490x330x400mm	5.0 Kg
D301.0313S	D301.0313C	9.0 Fr x 12 cm	05/30	318x235x128mm/490x330x400mm	5.0 Kg
D301.0314S	D301.0314C	10 Fr x 12 cm	05/30	318x235x128mm/490x330x400mm	5.0 Kg
D301.0325S	D301.0325C	11.5 Fr / 12 Fr x 13.5 cm	05/30	318x235x128mm/490x330x400mm	5.0 Kg
D301.0335S	D301.0335C	11.5 Fr / 12 Fr x 15 cm	05/30	318x235x128mm/490x330x400mm	5.0 Kg
D301.0345S	D301.0345C	11.5 Fr / 12 Fr x 16 cm	05/30	318x235x128mm/490x330x400mm	5.0 Kg
D301.0355S	D301.0355C	11.5 Fr / 12 Fr x 20 cm	05/30	318x235x128mm/490x330x400mm	5.0 Kg

Triple Lumen Catheter Subclavin & Curved Extensions

PRODUCT CODE					
Subclavin Extensions	Curved Extensions	Size	Standard Packing Box/Carton	Dimensions Box/Carton	Gross Wt.
T301.0325S	T301.0325C	11.5 Fr / 12 Fr x 13.5 cm	05/30	318x235x128mm/490x330x400mm	5.0 Kg
T301.0335S	T301.0335C	11.5 Fr / 12 Fr x 15 cm	05/30	318x235x128mm/490x330x400mm	5.0 Kg
T301.0345S	T301.0345C	11.5 Fr / 12 Fr x 16 cm	05/30	318x235x128mm/490x330x400mm	5.0 Kg
T301.0355S	T301.0355C	11.5 Fr / 12 Fr x 20 cm	05/30	318x235x128mm/490x330x400mm	5.0 Kg

Contents:

- One Polyurethane Radiopaque Cuffed Catheter with Injection Site Caps
- One Nitinol Core Guide Wire With Advancer
- Two Dilator
- One No. 11 Scalpel Blade
- One Introducer Needle
- One 5 ml Hollow Syringe
- Two Heparin Locks
- One IV Dressing
- One Surgical Drape
- One Mersilk Suture

SCEPTRE

REIN[®]
CATH

Single Lumen Hemodialysis Catheters

"Catheters and key components are developed based on research & are sourced from TAIWAN & PRC"

Clinically Tested

ISO 9001 : 2015 Certified

FDA Registered



Soft Tip

Reduces trauma to vessel wall upon introduction and in-situ.

Antimicrobial Heparin Cap

Specially-made adaptor helps to avoid infection.

Rotating Suture Wing

Rotating Suture Wing provides secure external anchoring.

Radiopaque Material

Sufficiently Radiopaque materials ensures the correct placement of catheter when visualized under X-Ray.

Biocompatible Polyurethane

Material tested for toxicity and hemocompatibility which becomes softer at body temperature.

Single Lumen Hemodialysis (Femoral) Catheter

Product Code	Size	Standard Packing Box/Carton	Dimensions Box/Carton
O301.0381S	7.0Fr (14G) x 13.0-14.0 cm	25/250	110x110x350mm/570x360x270mm
O301.0371S	8.0Fr (13G) x 13.0-14.0 cm	25/250	110x110x350mm/570x360x270mm
O301.0361S	9.0Fr (12G) x 13.0-14.0 cm	25/250	110x110x350mm/570x360x270mm

Note: Single Lumen Hemodialysis (Femoral) Catheter are also available with following components.

1. Dilator
2. Guidewire
3. Introducer Needle