



EndoView BR Series Single use Flexible Bronchoscope

Technical Specifications

Model		BR-M22	BR-M40A	BR-M50A	BR-M58A
Optical system	Field view	120°			
	Depth of field	3-50mm			
Insertion tube	Distal end outer diameter	2.3mm	4.1mm	5.1mm	6.1mm
	Working length	600mm			
Working Channel		/	2.0mm	2.8mm	3.4mm
Angulation	Bending	210° U/D			
Total Length		855mm			

Processor	
Function	Photo
	Video
	White Balance
	Freeze
Storage	Minimal 16G
Connection	DVI, HDMI
	VGA
	SDI, AV
Display	
Screen	LCD
	Size: 8 inch
Function	Photo
	Video
	White Balance
	Freeze
Battery	Type: Lithium
	Maximum 4 hours
Connection	HDMI, DVI


Shenzhen HugeMed
Medical Technical Development Co.,Ltd

Add: 416-1,516-1, Building 2, No. 1, Mawu Road, Baoan
Community, Yuanshan Street,Longgang District,
Shenzhen 518115, China
Tel: +86-755-2227 5866
+86-755-8226 7820
Fax: +86-755-2227 5833
Email: sale@hugemed.net

AN ENDOSCOPE
SOLUTION COMPANY

EndoView BR Series
Single use Flexible Bronchoscope


- Always Available
- Always Brand New
- No Costly Repair
- No Contamination




www.hugemed.net

EndoView BR Series Single use Flexible Bronchoscope

- Plug in and ready to go
- Multiple size to satisfy various clinical settings
- Navigate with ease




Standard **2 Million** Pixels,
High Quality Image



210° / 210° high deflection,
Ensure full navigation

- Slim: OD 3.2/ ID 1.2 mm
- Regular: OD 4.1/ ID 2.0 mm
- Regular: OD 5.1/ ID 2.8 mm
- Large: OD 6.1/ ID 3.4 mm



Compatible with conventional
accessories



Customized Two buttons:
WB / Photo / Video / Freeze

USB and HDMI connection,
easy data transfer