



Single-use Rhinolaryngoscope


Your ideal single-use rhinolaryngoscopy solution

HugeMed™ Single-use Rhinolaryngoscopes offer a safe, time-saving option for ENT examinations and airway assessments. Paired with the MS-8 featuring 15.6" FHD touchscreen, they help achieve high efficiency in imaging and clinical workflows across OPDs, operating rooms, and emergency settings.



Medical Image Processor

Image processing, display and recording



Specification



Model

MS-8

Dimension	386.4 (h) x 253.6 (w) x 38.5 (d) mm
Weight	2.7 kg
Display size	15.6 inches
Display resolution	1080P(1920*1080)px
Display type	TFT-LCD, touch screen
Function	Dual scope display, Brightness control, Contrast control, Using time calculator, etc.
Working time	Continuous working time ≤5 h
Memory capacity	Built-in SD card standard 64 GB

Single-use Rhinolaryngoscope

Great manoeuvrability and multiple sizes

Specification



Model	RH-M40H	RH-M52H	RH-M58H
Outer diameter	4.0 mm	5.0 mm	6.0 mm
Working channel	1.2 mm	2.2 mm	2.8 mm
Up/down bending angle	210°/210°	210°/210°	210°/210°
Field of view	120°	120°	120°
Working length	350 mm	350 mm	350 mm

Order Information

Product

Item Code



Single-use Rhinolaryngoscope	
RH-M40H	2018B-PA012
RH-M52H	2018B-PA011
RH-M58H	2018B-PA010
Medical Image Processor	
MS-8	2288B-PA001

Specifications, design and accessories are subject to change without any notice or obligation on the part of the manufacturer.

Shenzhen HugeMed Medical Technical Development Co., Ltd.

401, 501, Building 4, Haizhi Technology Park, Fortis, No.17, Bulan Road, Xialilang Community, Nanwan Street, Longgang District, Shenzhen, Guangdong, 518112, China

+86-755-22275833

+86-755-22275833

info@hugemed.net

hugemed.net

HugeMed