



IC-508
IV Stand

Light, agile
and easy to use

Features

Smooth Casters



The IV stand has five casters which roll smoothly with a small turning radius. Wheel-covers prevent the casters from slipping into narrow gaps.

Stable Base



A built-in step around the base can be supported with a foot to keep the IV stand from moving, a feature especially helpful during infusion bag replacement.

Steep Fall Prevention



A damper is provided inside the support pole to prevent sudden drops when the vertical stopper is loosened.

Quick and Easy Height Adjustment



Rotate the slide stopper to unlock, slide upper pole to desired height and re-rotate to lock in place.

Adaptable Grip Hold



The grip hold design has a characteristic shape and can be adjusted to 12 angles allowing adaption for any user.

Provision for Disinfection

The 1% sodium hypochlorite diluent can be used for chemical resistance and disinfection purpose such as during blood contamination.

Specifications

Model No.	IC-508
Dimensions	Base: 47.4 cm Total Height: 128.7 to 196.2 cm
Product Weight	4.3 kg
Max Load Weight	Hooks: 1.5 kg (per hook) 2.5 kg (both hooks combined) Sliding pipe (upper pole): 3.5 kg Pole (upper + lower pole): 12.5 kg
Material	Pole: Stainless Steel Base, Grip Hold, Vertical Stopper: Polypropylene

