

HugeMed™

Airway Management Workstation

One trolley, equipped with all the airway management visual tools. Base on the universal compatibility, the MS-8 can be connected to the VLRM/VLDM/VL3S/VL4REX, as well as all the single-use bronchoscope to achieve one solution for airway management.

Video Laryngoscope



Single-use Bronchoscope



Image Processor MS-8



HugeMed

Shenzhen HugeMed Medical Technical Development Co., Ltd.
415, 416-1, 516-1, Building 2, No.1, Mawu Road, Baoan Community, Yuanshan Street,
Longgang District, Shenzhen, Guangdong, 518115, China
☎ +86-755-22275833 ☎ +86-755-22275833 ✉ info@hugemed.net

For more information
please visit our website:
www.hugemed.net


CE

P/N: ENG-EU-CATA-IP-MS-8-202309V1.0

HugeMed

MS-8 15.6 inches touch screen Image Processor

The all-around upgrade digital video processor for HugeMed™ adopts a new operating system, which functions as a central system to facilitate easier reading and operating in various endoscopy applications. Especially important for occasions where mobility is needed in hospitals.




Dual Scopes Display

The dual scopes display interface allows the views of the flexible type scope and blade type scope simultaneously, which enhances the efficiency during the intubation procedures.



Diagnostic Report


The built-in diagnostic report generation system enables the output of diagnostic reports with one click to increase the flow efficiency of endoscopy.



IT Solution

With DICOM protocol, the MS-8 recorded image and video can be transferred to hospital PACS system. Effectively improving the efficiency of inter-departmental collaboration.


Multipurpose Image processor



Suitable for a variety of clinical allocation scenario.

Compatibility

- Gold fingers and Type-C connector for scope
- More efficient data transfer capabilities
- Plug and unplug more smoothly
- Effectively reduce plug damage



Gold fingers connector



TYPE-C



More Connection abilities

- USB3.0 for faster data transfer
- Wired or wireless networks ability
- DVI/HDMI/SDI image output ability



USB



HDMI



DVI



Wifi



LAN



Flexible Mobility

Besides being used on the table, the MS-8 can be equipped with the trolley for different occurrences such as OR, ICU, OPD, etc.

